

Product Datasheet

Bovine Insulin-like Growth Factors 2 ELISA Kit (orb1214712)

Description	Bovine Insulin-like Growth Factors 2, IGF-2 ELISA Kit
Reactivity	Bovine
Range	15-3000ng/ml
Tested Applications	ELISA
Target	IGF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,066
Uniprot ID	P07455
Solubility (25°C)	3200ng/ml
Sensitivity	7.19ng/ml
Entrez	281239
Expiration Date	Please enquire.