

Product Datasheet

Bovine Peroxisome Proliferation Activated Receptor Gamma ELISA Kit (orb1214713)

Description	Bovine Peroxisome Proliferation Activated Receptor Gamma, PPARG ELISA Kit
Reactivity	Bovine
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	PPARG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,579
Uniprot ID	018971
Solubility (25°C)	24ng/ml
Sensitivity	0.016ng/ml
Entrez	281993
Expiration Date	Please enquire.