

Product Datasheet

Bovine Interleukin-2 Receptor Subunit Alpha ELISA Kit (orb1214728)

Description	Bovine Interleukin-2 Receptor Subunit Alpha, IL2RA ELISA Kit
Reactivity	Bovine
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	IL2RA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,239
Uniprot ID	P12342
Solubility (25°C)	48ng/mL
Sensitivity	0.048ng/mL
Entrez	281861
Expiration Date	Please enquire.