

## Product Datasheet

### Bovine Appetite-Regulating Hormone ELISA Kit (orb1214743)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Bovine Appetite-Regulating Hormone, GHRL ELISA Kit |
| <b>Reactivity</b>          | Bovine                                             |
| <b>Range</b>               | 0.3-90ng/L                                         |
| <b>Tested Applications</b> | ELISA                                              |
| <b>Target</b>              | GHRL                                               |
| <b>Storage</b>             | 2-8°C                                              |
| <b>Note</b>                | For research use only                              |
| <b>Application notes</b>   | Sandwich                                           |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates             |
| <b>MW</b>                  | 12,793                                             |
| <b>Uniprot ID</b>          | <b>Q9BDJ6</b>                                      |
| <b>Solubility (25°C)</b>   | 96ng/L                                             |
| <b>Sensitivity</b>         | 0.16ng/L                                           |
| <b>Entrez</b>              | <b>281192</b>                                      |
| <b>Expiration Date</b>     | Please enquire.                                    |