

Product Datasheet

Bovine Kappa-Casein ELISA Kit (orb1214767)

Description	Bovine Kappa-Casein, CSN3 ELISA Kit
Reactivity	Bovine
Range	7-1500ug/ml
Tested Applications	ELISA
Target	CSN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21269
Uniprot ID	P02668
Solubility (25°C)	1600µg/ml
Sensitivity	3.33µg/ml
Entrez	281728
Expiration Date	Please enquire.