

Product Datasheet

Bovine T-cell Surface Glycoprotein Cd8 Alpha Chain ELISA Kit (orb1214769)

Description	Bovine T-cell Surface Glycoprotein Cd8 Alpha Chain, CD8A ELISA Kit
Reactivity	Bovine
Range	2-600ng/ml
Tested Applications	ELISA
Target	CD8A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,417
Uniprot ID	P31783
Solubility (25°C)	800ng/ml
Sensitivity	0.91ng/ml
Entrez	281060
Expiration Date	Please enquire.