

Product Datasheet

Bovine Interleukin 18 ELISA Kit (orb1214772)

Description	Bovine Interleukin 18, IL-18 ELISA Kit
Reactivity	Bovine
Range	20-4500ng/L
Tested Applications	ELISA
Target	IL18
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,347
Uniprot ID	Q9TU73
Solubility (25°C)	4800ng/L
Sensitivity	9.99ng/L
Entrez	281249
Expiration Date	Please enquire.