

Product Datasheet

Bovine Gastrin-releasing Peptide ELISA Kit (orb1214777)

Description	Bovine Gastrin-releasing Peptide, GRP ELISA Kit
Reactivity	Bovine
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	GRP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,934
Uniprot ID	Q863C3
Solubility (25°C)	96ng/ml
Sensitivity	0.16ng/ml
Entrez	615323
Expiration Date	Please enquire.