

Product Datasheet

Bovine Interleukin 2 ELISA Kit (orb1214781)

Description	Bovine Interleukin 2, IL-2 ELISA Kit
Reactivity	Bovine
Range	1-400ng/L
Tested Applications	ELISA
Target	IL2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,628
Uniprot ID	P05016
Solubility (25°C)	480ng/L
Sensitivity	0.47ng/L
Entrez	280822
Expiration Date	Please enquire.