

Product Datasheet

Bovine Deoxyribonuclease 1 ELISA Kit (orb1214789)

Description	Bovine Deoxyribonuclease 1, DNASE-1 ELISA Kit
Reactivity	Bovine
Range	7-1500ng/L
Tested Applications	ELISA
Target	DNASE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,346
Uniprot ID	P00639
Solubility (25°C)	1600ng/L
Sensitivity	3.91ng/L
Entrez	282217
Expiration Date	Please enquire.