

Product Datasheet

Bovine Gonadotropin-releasing Hormone Receptor ELISA Kit (orb1214791)

Description	Bovine Gonadotropin-releasing Hormone Receptor, GNRHR ELISA Kit
Reactivity	Bovine
Range	30-9000pg/ml
Tested Applications	ELISA
Target	GNRHR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,741
Uniprot ID	P32236
Solubility (25°C)	9600pg/ml
Sensitivity	15pg/ml
Entrez	281798
Expiration Date	Please enquire.