

Product Datasheet

Bovine Atrial Natriuretic Peptide Receptor 2 ELISA Kit (orb1214796)

Description	Bovine Atrial Natriuretic Peptide Receptor 2, NPR2 ELISA Kit
Reactivity	Bovine
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	NPR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	117,146
Uniprot ID	P46197
Solubility (25°C)	16ng/ml
Sensitivity	0.023ng/ml
Entrez	281357
Expiration Date	Please enquire.