

Product Datasheet

Bovine Endothelin 1 ELISA Kit (orb1214797)

Description	Bovine Endothelin 1, ET-1 ELISA Kit
Reactivity	Bovine
Range	2-600ng/L
Tested Applications	ELISA
Target	EDN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,977
Uniprot ID	P17322
Solubility (25°C)	640ng/L
Sensitivity	1.19ng/L
Entrez	281137
Expiration Date	Please enquire.