

Product Datasheet

Bovine Apolipoprotein H ELISA Kit (orb1214803)

Description	Bovine Apolipoprotein H, APO-H ELISA Kit
Reactivity	Bovine
Range	0.2-70ug/ml
Tested Applications	ELISA
Target	APOH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,252
Uniprot ID	P17690
Solubility (25°C)	80µg/ml
Sensitivity	0.11µg/ml
Entrez	281006
Expiration Date	Please enquire.