

Product Datasheet

Bovine Interleukin 4 ELISA Kit (orb1214812)

Description	Bovine Interleukin 4, IL-4 ELISA Kit
Reactivity	Bovine
Range	7-1500ng/L
Tested Applications	ELISA
Target	IL4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,119
Uniprot ID	P30367
Solubility (25°C)	1600ng/L
Sensitivity	3.34ng/L
Entrez	280824
Expiration Date	Please enquire.