

Product Datasheet

Bovine Interferon Gamma ELISA Kit (orb1214816)

Description	Bovine Interferon Gamma, IFN-G ELISA Kit
Reactivity	Bovine
Range	5-2000pg/ml
Tested Applications	ELISA
Target	IFNG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,393
Uniprot ID	P07353
Solubility (25°C)	2400pg/ml
Sensitivity	2.35pg/ml
Entrez	281237
Expiration Date	Please enquire.