

## Product Datasheet

### CD3 antibody (orb1215079)

|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD3                                                                  |
| <b>Species/Host</b>        | Mouse                                                                                             |
| <b>Reactivity</b>          | Human                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                      |
| <b>Tested Applications</b> | FC                                                                                                |
| <b>Storage</b>             | 2-8°C                                                                                             |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | Mouse IgG2a, kappa                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">P20963</a> , <a href="#">P09693</a> , <a href="#">P07766</a> , <a href="#">A8MVP6</a> |
| <b>Entrez</b>              | <a href="#">915</a>                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |