

Product Datasheet

CD45RO antibody (Biotin) (orb1215475)

Description	Mouse monoclonal antibody against CD45RO conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IF, IHC-Fr, IHC-P, IP, WB
Storage	2-8°C protected from light
Note	For research use only
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P08575
Entrez	5788
Expiration Date	12 months from date of receipt.