

Product Datasheet

CD45RA antibody (APC/Cy7) (orb1215509)

Description	Mouse monoclonal antibody against CD45RA conjugated to APC/Cy7
Species/Host	Mouse
Reactivity	Human
Conjugation	APC/Cy7
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P08575
Entrez	5788
Expiration Date	12 months from date of receipt.