

Product Datasheet

CD3 antibody (APC) (orb1215632)

Description	Mouse monoclonal antibody against CD3 conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P20963 , P09693 , P07766 , A8MVP6
Entrez	915
Expiration Date	12 months from date of receipt.