


## Product Datasheet

### MPXV A29L Protein (orb1215650)

|                            |                                                                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Recombinant MPXV A29L Protein is produced by E.coli expression system and the target gene encoding S...                                          |
| <b>Reactivity</b>          | Virus                                                                                                                                            |
| <b>Tested Applications</b> | ELISA, MS, SDS-PAGE, WB                                                                                                                          |
| <b>Preservatives</b>       | Liquid in a 0.2 µm filtered solution of PBS, 300 mM NaCl, pH 7.4.                                                                                |
| <b>Form/Appearance</b>     | Liquid in a 0.2 µm filtered solution of PBS, 300 mM NaCl, pH 7.4.                                                                                |
| <b>Storage</b>             | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                            |
| <b>Application notes</b>   | Please minimize freeze-thaw cycles.                                                                                                              |
| <b>Purity</b>              | Greater than 95% as determined by reducing SDS-PAGE.                                                                                             |
| <b>Source</b>              | E.coli                                                                                                                                           |
| <b>Uniprot ID</b>          | <b>P48865</b>                                                                                                                                    |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                   |



Greater than  
95% as  
determined  
by reduci...