

Product Datasheet

Guinea pig TNF alpha protein (orb1215677)

Description	The Guinea Pig TNF alpha Biotinylated yeast-derived recombinant protein is not tagged and is natural...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SASQNDNDKP VAHVVANQQA EEELQWLSKR ANALLANGMG LSDNQLVVP DGLYLIYSQV LFKGQGCPY LLLTHTVSRL AVSYPEKVN LSAIKSPCK ETPEGAERK WYEPIYLG V FQLQKGDRL AEVNLPQYLD FADSGQIYFG VIAL (154)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	100135630
Expiration Date	6 months from date of receipt.