

Product Datasheet

Human IL-2 protein (orb1216191)

Description	The Human IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTKKT QLQLEHLLLD LQMILNGINN YKNPKLTRML TFKFYMPKKA TELKHLQCLE EELKPLEEVL NLAQSKNFHL RPRDLISNIN VIVLELKGSE TTFMCEYADE TATIVEFLNR WITFCQSIIS TLT (133)
Purity	98%
Source	Yeast
MW	15.4 kDa
Entrez	3558
Expiration Date	6 months from date of receipt.