

Product Datasheet

Canine CXCL10 protein (orb1216195)

Description	The Canine CXCL10 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CXCL10
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IPLSRTIRCT CIKISDQPVN LRSLEKIEMI PASPSCPHVE IIATMKKSGE KRCLNPESKT IKSLVKAISK KRSRRSP (77)
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	478432
Expiration Date	6 months from date of receipt.