

Product Datasheet

Mouse B7-H2 protein (orb1216197)

Description	The Mouse B7-H2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	B7-H2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ETEVGAMVGS NVVLSCIDPH RRHFNLSGLY VYWQIENPEV SVTTYLPYKS PGINVDSSYK NRGHLSLDSM KQGNFSLYLK NVTPQDTQEF TCRVFMNTAT ELVKILEEVV RLRVAANFST PVISTSDSSN PGQERTYTCM SKNGYPEPNL YWINTTDNSL IDTALQNNTV YLNKLGLYDV ISTLRLPWTS HGDVLCCVEN VALHQNITSI SQAESFTGNN TKNPQETHNN ELKVLVPVLA VLAAAAFVSF IYRRTRPHR SYTGPKTVQL ELTDHA (276) [His185Arg]
Purity	98%
Source	Yeast
MW	31.0 kDa
Entrez	50723
Expiration Date	6 months from date of receipt.