

Product Datasheet

Mouse IL-13 protein (orb1216200)

Description	The Mouse IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPVPRSVSLP LTLKELIEEL SNITQDQTPL CNGSMVWSVD LAAGGFCVAL DSLTNISNCN AIYRTQRILH GLCNRKAPTT VSSLPDTKIE VAHFITKLLS YTKQLFRHGP F (111)
Purity	98%
Source	Yeast
MW	12.2 kDa
Entrez	16163
Expiration Date	6 months from date of receipt.