

## Product Datasheet

### Mouse M-CSF protein (orb1216220)

|                         |                                                                                                                                                                                                                                                                                          |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | The Mouse M-CSF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...                                                                                                                                                                                  |
| <b>Target</b>           | M-CSF                                                                                                                                                                                                                                                                                    |
| <b>Form/Appearance</b>  | Lyophilized                                                                                                                                                                                                                                                                              |
| <b>Storage</b>          | -20°C                                                                                                                                                                                                                                                                                    |
| <b>Note</b>             | For research use only                                                                                                                                                                                                                                                                    |
| <b>Protein Sequence</b> | AKEVSEHCSH MIGNGHLKVL QQLIDSQMET<br>SCQIAFEFVD QEQLDDPVCY LKKAFFLVQD<br>IIDETMRFKD NTPNANATER LQELSNNLNS<br>CFTKDYEEQN KACVRTFHET PLQLEKIKN<br>FFNETKNLLE KDOWNIFTKNC NNSFAKCSSR<br>DVVTKPDCNC LYPKATPSSD PASASPHQPP<br>APSMAPLAGL AWDDSQRTEG SLLPSELPL<br>RIEDPGSAKQ RPPRSTCQTL E (231) |
| <b>Purity</b>           | 98%                                                                                                                                                                                                                                                                                      |
| <b>Source</b>           | Yeast                                                                                                                                                                                                                                                                                    |
| <b>MW</b>               | 26.0 kDa                                                                                                                                                                                                                                                                                 |
| <b>Entrez</b>           | <a href="#">12977</a>                                                                                                                                                                                                                                                                    |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                                                                                                                                                                                                           |