

Product Datasheet

Human TNF alpha protein (orb1216226)

Description	The Human TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VRSSSRTPSD KPVAVVAVNP QAEGQLQWLN RRANALLANG VELRDNQLVV PSEGLYLIYS QVLFKGQGCP STHVLLHTI SRIAVSYQTK VNLLSAIKSP CQRETPEGAE AKPWYEPIYL GGVFQLEKGD RLSAENRPDY LDFAESGQVY FGIIAL (156) [Ile135/136Deletion]
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	7124
Expiration Date	6 months from date of receipt.