

Product Datasheet

Canine IL-8 protein (orb1216335)

Description	The Canine IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VSELRCQCIKTHSTPFHPKYIKELRVIDSGPHCENSEIIVKLFNGNEVCLDPKEKW VQKVVFIFLKKAQKQDP
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	403850
Expiration Date	6 months from date of receipt.