

## Product Datasheet

### Mouse IL-22 ELISA Kit (orb1216560)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The Mouse IL-22 ELISA contains capture antibody, standard, and detection antibody for development of... |
| <b>Reactivity</b>          | Mouse                                                                                                   |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Target</b>              | IL-22                                                                                                   |
| <b>Storage</b>             | 4 degrees C.                                                                                            |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |