

Product Datasheet

IL-8 antibody (orb1216667)

Description	The Ferret IL-8 (CXCL8) polyclonal antibody is unlabeled and has been qualified for use in ELISA and...
Species/Host	Rabbit
Reactivity	Ferret
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Ferret IL-8
Target	IL-8
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.