

Product Datasheet

Fas Ligand antibody (orb1216901)

Description	The Swine Fas Ligand polyclonal antibody is unlabeled and has been qualified for use in ELISA and We...
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Swine Fas Ligand
Target	Fas Ligand
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.