

Product Datasheet

IFN alpha antibody (orb1216959)

Description	The Chicken IFN alpha polyclonal antibody is unlabeled and has been qualified for use in ELISA and W...
Species/Host	Goat
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Chicken IFN alpha
Target	IFN alpha
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.