

## Product Datasheet

### BLZ945 (orb1223530)

|                        |                                                       |
|------------------------|-------------------------------------------------------|
| <b>Description</b>     | BLZ945                                                |
| <b>Storage</b>         | Storage temperature: -20°C. Stability: $\geq 2$ years |
| <b>Note</b>            | For research use only                                 |
| <b>Purity</b>          | > 98% (HPLC)                                          |
| <b>CAS Number</b>      | 953769-46-5                                           |
| <b>Expiration Date</b> | 12 months from date of receipt.                       |