

## Product Datasheet

### KX2-391 (orb1223723)

|                        |                                                       |
|------------------------|-------------------------------------------------------|
| <b>Description</b>     | KX2-391                                               |
| <b>Storage</b>         | Storage temperature: -20°C. Stability: $\geq 2$ years |
| <b>Note</b>            | For research use only                                 |
| <b>Purity</b>          | > 98%(HPLC)                                           |
| <b>CAS Number</b>      | 897016-82-9                                           |
| <b>Expiration Date</b> | 12 months from date of receipt.                       |