

## Product Datasheet

### Milrinone (orb1224003)

|                        |                                                  |
|------------------------|--------------------------------------------------|
| <b>Description</b>     | Milrinone                                        |
| <b>Storage</b>         | Storage temperature: -20°C. Stability: ≥ 2 years |
| <b>Note</b>            | For research use only                            |
| <b>Purity</b>          | > 98%(HPLC)                                      |
| <b>CAS Number</b>      | 78415-72-2                                       |
| <b>Expiration Date</b> | 12 months from date of receipt.                  |