



Product Datasheet

Dnajc5 antibody (orb12269)

Description	Goat polyclonal antibody to Dnajc5
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Dnajc5 (aa177-191)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 20kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 22.1kDa according to NP_058055.1). Recommended concentration: 0.5-1.5µg/ml. Preliminary testing was unsuccessful on Mouse and Rat Brain for this particular batch.. IHC: Paraffin embedded Human Skeletal Muscle. Recommended concentration: 5µg/ml.
Protein Sequence	SATETTQLTADSHPS
Clonality	Polyclonal
MW	22.1
Entrez	80331
Dilution Range	ELISA: 1:64000, WB: 0.5-1.5 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb12269 (0.5µg/ml)
staining of Human
C...



Immunohistochemical
staining of Human
Sk...