

Product Datasheet

SARS-CoV-2 trimeric soluble full-length Spike protein, Omicron variant, S protein (orb1227191)

Description	Recombinant SARS-CoV-2 trimeric soluble full-length Spike protein, Omicron variant, S protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	8xHis at the C-terminus
Note	For research use only
Protein Sequence	16-1211
Purity	> 95%
Source	HEK293 cells
MW	137
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.