

## Product Datasheet

### HeLa Cell Slide [Out of stock] (orb1227248)

<b>Description</b>	HeLa Cell Slide [Out of stock]
<b>Tested Applications</b>	ICC
<b>Storage</b>	Hela cell slides can be stored for up to 6 months at 4°C under desiccate conditions.
<b>Note</b>	For research use only
<b>Application notes</b>	Hela cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Before applying immunocytochemistry, 0.1% Triton X-100 pretreatment for 1 minute at room temperature is highly suggested. Optimal conditions should be determined by the end user.
<b>Dilution Range</b>	Hela cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Before applying immunocytochemistry, 0.1% Triton X-100 pretreatment for 1 minute at room temperature is highly suggested. Optimal conditions should be determined by the end user.
<b>Expiration Date</b>	12 months from date of receipt.