

Product Datasheet

SARS-CoV-2 Omicron BA.1 Variant Pseudovirus (orb1227402)

Description	SARS-CoV-2 Omicron BA.1 Variant Pseudovirus
Tested Applications	NeA
Target	S
Preservatives	DMEM/10% FBS
Form/Appearance	Liquid
Storage	The pseudovirus is stable in a liquid state at -70° for 12 months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Protein Sequence	Omicron Variant Spike BA.1 Met 1 - Cys 1255 (A67V, HV69-70del, T95I, G142D, VYY143-145del, N211del, L212I, ins214EPE, G339D, S371L, S373P, S375F, K417N, N440K, G446S, S477N, T478K, E484A, Q493R, G496S, Q498R, N501Y, Y505H, T547K, D614G, H655Y, N679K, P681H, N764K, D796Y, N856K, Q954H, N969K, L981F)
Source	HEK293 cells
Uniprot ID	P0DTC1
NCBI	QHD43416
Expiration Date	12 months from date of receipt.