

Product Datasheet

Adipokines: Adiponectin and CTRPs 1-7 Detection Set (orb1227481)

Description	Adipokines: Adiponectin and CTRPs 1-7 Detection Set
Reactivity	Human
Tested Applications	IF, IHC, WB
Immunogen	Rabbit polyclonal antibodies were raised against peptides corresponding to amino acid sequences from each of the corresponding proteins.
Preservatives	PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	Antibody 1 mg/mL
Storage	Stable at 4°C for three months, store at -20°C for up to one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	Please enquire.