

Product Datasheet

Clevudine (orb1227508)

Description	Clevudine
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4°C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	C ₁₀ H ₁₃ FN ₂ O ₅
Purity	greater than or equal to 98% (HPLC)
MW	260.2
CAS Number	163252-36-6
SMILES	<chem>O[C@H]1[C@H](CO)O[C@]([H])(N2C=C(C)C(NC2=O)=O)[C@@H]1F</chem>
Solubility (25°C)	Soluble in water (30mg/ml) or DMSO (20mg/ml).
Expiration Date	12 months from date of receipt.

