

Product Datasheet

Favipiravir (orb1227532)

Description	Favipiravir
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4C. Long Term Storage: -2°C. Handling Advice: Keep cool and dry. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	C ₅ H ₄ FN ₃ O ₂
Purity	greater than or equal to 98%
MW	157.1
CAS Number	259793-96-9
SMILES	FC1=CNC(C(C(N)=O)=N1)=O
Solubility (25°C)	Soluble in water (30mg/ml) or DMSO (30mg/ml).
Expiration Date	12 months from date of receipt.

