

Product Datasheet

Ebselen (orb1227536)

Description	Ebselen
Form/Appearance	Lyophilized
Storage	Short Term Storage: +4°C. Long Term Storage: -2°C. Handling Advice: Protect from light. Use/Stability: Stable for at least 2 years after receipt when stored at -2°C.
Note	For research use only
Formula	C ₁₃ H ₉ NOSe
Purity	greater than or equal to 98% (HPLC)
MW	274.2
CAS Number	60940-34-3
SMILES	<chem>O=C1N([Se]C2=CC=CC=C2)C1=CC=CC=C1</chem>
Solubility (25°C)	Soluble in DMSO (25mg/ml), DMF (25mg/ml) or water (2mg/ml).
Expiration Date	12 months from date of receipt.

