


Product Datasheet

PRDM14 antibody (orb12276)

Description	Goat polyclonal antibody to PRDM14
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRDM14 (aa533-544)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in lysates of cell line A549 (calculated MW of 64.1kDa according to NP_078780.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	DSHVRRSHKEDD
Clonality	Polyclonal
MW	64.1
Entrez	63978
Dilution Range	ELISA: 1:2000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate using PRDM14 antibody (orb12276).