

Product Datasheet

SARS-CoV-2 (COVID-19) S1 Recombinant Protein (orb1227613)

Description SARS-CoV-2 (COVID-19) S1 Recombinant Protein

Tested Applications WB

Target S

Preservatives Lyophilized from 0.22 um filtered solution in Tris with Glycine, Arginine and NaCl, pH 7.5. Normally trehalose is added as protectant before lyophilization.

Form/Appearance Lyophilized

Storage Lyophilized protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Tag Fc Tag

Note For research use only

Application notes This recombinant protein can be used for WB. For research use only.

Protein Sequence Val 16 - Arg 685

Purity > 95% as determined by SDS-PAGE.


Source HEK293 cells

MW 101.5 kDa

NCBI [QHD43416.1](#)

Dilution Range This recombinant protein can be used for WB. For research use only.

Expiration Date 6 months from date of receipt.



SARS-CoV-2
(COVID-19)
S1
recombinant
pro...