


Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Recombinant Protein (orb1227614)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Recombinant Protein
Tested Applications	WB
Target	NS1
Preservatives	Lyophilized from 0.22 um filtered solution in PBS, pH 7.4. Normally trehalose is added as protectant before lyophilization.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Met 1 - Ala 419
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	49.4 kDa
NCBI	QH062115.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



2019-nCoV
nucleocapsid
recombinant
prote...