

Product Datasheet

SARS-CoV-2 (COVID-19) Envelope Recombinant Protein (orb1227615)

Description	SARS-CoV-2 (COVID-19) Envelope Recombinant Protein
Target	E
Preservatives	Lyophilized from 0.22 um filtered solution in PBS, pH 7.4. Normally trehalose is added as protectant before lyophilization.
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	GST Tag, His Tag
Note	For research use only
Protein Sequence	Met 1 -Val 75
Source	E. coli
MW	9.6 kDa
NCBI	QHU79206.1
Expiration Date	6 months from date of receipt.