


Product Datasheet

VEGF 120 Recombinant Protein (HEK293 Expressed) (orb1227838)

Description	VEGF 120 Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	Vegfa
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Ala 27 - Arg 146
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	17.2 kDa
Uniprot ID	Q00731
NCBI	AAB22254.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Mouse
VEGF120,
Tag Free
on SDS-
PAGE
unde...