


Product Datasheet

UCHL1 Recombinant Protein (orb1227842)

Description	UCHL1 Recombinant Protein
Tested Applications	WB
Target	UCHL1
Preservatives	20mM Tris, 500mM NaCl, 1mM DTT, pH 8.0
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -71°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Met 1 - Ala 223
Purity	> 92% as determined by SDS-PAGE.
Source	E. coli
MW	25.7 kDa
Uniprot ID	P09936
NCBI	NP_004172
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of
rh UCHL1
was
determined
by...